**EN INSTRUCTIONS OF USE**

SoftFil® Precision

DESCRIPTION

SoftFil® Precision are sterile and single use injection kits containing:
- 1 stainless steel pre-hole needle in its protective cap.
- 1 stainless steel, graduated blunt-tip needle (micro-cannula) in its protective cap.
SoftFil® micro-cannulas are available in various lengths and gauge.

INDICATIONS

SoftFil® micro-cannulas are intended for hypodermal injection of body fluids or tissues, drugs or medical devices approved for human use.

CONTRAINDICATIONS

The use of micro-cannulas is contraindicated for vascular, medullar or nerve injections.

PRECAUTIONS

This device can be used only by or on the order of a physician and/or authorized healthcare professionals in accordance with local legislation and also trained in injection techniques. Before use, users shall inspect the injection site to make sure that there is no infection, inflammation, oedema, skin lesion, scar, birthmark, bony prominences or lipothrophy. The micro-cannula is chosen based on two main parameters which determine how flexible the micro-cannula is:
- the length, according to the extent of the area to be treated,
- the outside caliber (outside diameter), indicated in gauge (G) and in millimeters (mm), to be adapted to the viscosity of the product and to the injection plane.

The technique, method and products used depend on experience and are left to the physician and/or authorized healthcare professional's discretion in order to ensure the best clinical results. In case of significant bleeding at entry point, it is recommended to make a new entry point. Users should perform under appropriate aseptic conditions and handle needles with precaution to avoid inadvertent sharps injury and getting infectious disease.

Micro-cannulas are fragile devices, to handle with caution. Any deflection, crushing or damage of the tube may lead to its breaking during its use. Avoid any sudden movement or force in case of resistance. In case of breaking during the use, the patient should be addressed as soon as possible to an authorized surgeon who will take the appropriate measures. Do not use if the product or the pouch appears to be damaged. For single patient use only - Contents supplied sterile

Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in injury or illness. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause infection or cross infection including, but not limited to, the transmission of infectious diseases(s). Any device which has been used or partially used must be destroyed in accordance with the local regulations concerning the disposal of hospital waste.

USE

- 1/ Inspect the packaging to confirm its integrity before opening.
- 2/ Prepare the skin area to be injected: make-up removal, cleaning, rinsing, drying, antiseptic to prevent a biofilm.
- 3/ Open the packaging.
- 4/ To fix the SoftFil® micro-cannula on the syringe, hold the micro-cannula by its plastic cover firmly and:

- in case of syringes with fixed female thread: screw it towards the syringe tip.
- in case of syringes with mobile luer lock: push it strongly towards the syringe tip while twisting it clockwise. For optimal screwing, use a clamp or if necessary.
- 5/ Create pre-hole at the identified entry point, with the pre-hole needle included in the sterile injection kit.
- 6/ Through this entry hole, insert with one hand the micro-cannula previously securely connected and locked to the syringe into the desired injection level. To do so, pinch the skin between the thumb and the forefinger of the opposite hand.

Rely on the mark located on the hub indicating the orientation of the side port and the graduations every centimeter along the tube for determining the depth of the injection to favor a homogeneous treatment.

Gently proceed with the injection, applying the least amount of pressure.

7/ To remove the SoftFil® micro-cannula from the syringe, hold the syringe and turn the mobile female luer lock counterclockwise firmly, then pull the micro-cannula backwards gently.

Dispose of micro-cannulas and pre-hole needles in proper sharps container.

SIDE EFFECTS

Few complications as pain, redness, bruising, hematoma, swelling or oedema, may result from the injection technique itself. Furthermore, repeated injections can cause tissue scarring.

VALIDITY / STORAGE CONDITIONS

Do not use after the expiry date. The product must be stored in its original packaging, in a dry place and in the shelter of light.

SYMBOLS USED ON PACKAGING

REF	Product reference	Use until	LOT	Batch number	Use until	LOT	Batch number	Use until	LOT	Batch number
Do not reuse	Do not resterilize	Do not use if package is opened or damaged								
STERILE	Sterile product. Method of sterilization – Ethylene oxide	Dimension								
See the Instructions for use leaflet										
Caution:	Federal Law (USA) restricts this device to sale by or on the order of a physician									
Manufacturer	Keep away from moisture	Keep away from sunlight								
Micro-cannula	Pre-hole needle	Latex free								
0459 CE Mark	obtained in compliance with Directive 93/42/EEC related to medical devices.									
0459 is the number of the Notified Body.	CE Mark obtained in 2011.									

COLOUR CODES ON HUBS AND LABELS

Light green = 14G	White = 16G	Pink = 18G	Yellow = 20G	Black = 22G
Dark blue = 23G	Medium purple = 24G	Orange = 25G	Brown = 26G	
Yellow = 30G				

All serious adverse events associated with use of this device should be reported to the manufacturer. All returned products must first be decontaminated, and accompanied by a decontamination certificate. If no certificate is provided, the claim will be rejected.

Soft Medical Aesthetics 55 boulevard Periere – 75017 Paris – France Tel. +33 170 629 009 Email : info@softfil.com www.softfil.com	Soft Medical Aesthetics 55 boulevard Periere – 75017 Paris – France Tel. +33 170 629 009 Email : info@softfil.com www.softfil.com
--	--

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REF	Référence catalogue	Date limite d'utilisation	LOT	Code de lot	Date de fabrication
Do not reuse	No resterilize	Do not use if package is opened or damaged			
STERILE	Stérilisé avec de l'oxyde d'éthylène	Dimensions			
See the Instructions for use leaflet					
Caution:	Federal Law (USA) restricts this device to sale by or on the order of a physician				
Manufacturer	Keep away from moisture	Keep away from sunlight			
Micro-cannula	Pre-hole needle	Latex free			
0459 CE Mark	obtained in compliance with Directive 93/42/CEE related to medical devices.				
0459 is the number of the Notified Body.	CE Mark obtained in 2011.				

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