Soft Fil® Precision

**DESCRIPTION**

SoftFil® Precision are sterile and single use injection kits containing:

- 1 stainless steel pre-hole needle in its protective cap.
  - 1 stainless steel, graduated blunt-tip needle (micro-cannula) in its protective cap.
- SoftFil® micro-cannulas are available in various lengths and gauge.

**INDICATIONS**

SoftFil® micro-cannulas are intended for hypodermal injection of body fluids or tissues, drugs or medical devices approved for human use.

**CONTRAINDICATIONS**

The use of micro-cannulas is contraindicated for vascular, medullar or nerve injections.

**PRECAUTIONS**

This device can be used only by or on the order of a physician and/or authorized healthcare professionals in accordance with local legislation and also trained in injection techniques.

Before use, users shall inspect the injection site to make sure that there is no infection, inflammation, oedema, skin lesion, scar, birthmark, bony prominences or lipohypertrophy.

The micro cannula is chosen based on two main parameters which determine how flexible the micro-cannula is:

- the length, according to the extent of the area to be treated.
- the outside caliber (outside diameter), indicated in gauge (G) and in millimeters (mm), to be adapted to the viscosity of the product and to the injection plane.

The technique, method and products used depend on experience and are left to the physician and/or authorized healthcare professional's discretion in order to ensure the best clinical results.

In case of significant bleeding at entry point, it is recommended to make a new entry point.

Users should prepare under appropriate aseptic conditions and handle needles with precaution to avoid inadvertent sharps injury and getting infectious disease.

Micro-cannulas are fragile devices, to handle with caution. Any deflection, crushing or damage of the tube may lead to its breaking during its use. Avoid any sudden movement or force in case of resistance.

In case of breaking during the use, the patient should be addressed as soon as possible to an authorized surgeon who will take the appropriate measures.

Do not use if the product or the pouch appears to be damaged.

For single patient use only - Contents supplied sterile with ethylene oxide.

Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and lead to device failure, which in turn may result in injury or illness.

Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause infection or cross infection including, but not limited to, the transmission of infectious disease(s).

Any device which has been used or partially used must be destroyed in accordance with the local regulations concerning the disposal of hospital waste.

**USE**

- 1/ Inspect the packaging to confirm its integrity before opening.
- 2/ Prepare the skin area to be injected: make-up removal, cleaning, rinsing, drying, antiseptic to prevent a biofilm.
- 3/ Open the packaging.
- 4/ To fix the SoftFil® micro-cannula on the syringe, hold the micro-cannula by its plastic cover firmly and:
  - in case of syringes with fixed female thread: screw it towards the syringe tip.
  - in case of syringes with mobile female luer lock: push it strongly towards the syringe tip while twisting it clockwise. For optimal screwing, use a clamp if necessary.
- 5/ Create pre-hole at the identified entry point, with the SoftFil® Precision pre-hole needle included in the sterile injection kit.
- 6/ Through this entry hole, insert the micro-cannula previously securely connected and locked to the syringe with one hand by sliding the micro-cannula along the SoftFil® EasyGuide pre-hole needle into the desired injection level. Then, remove the pre-hole needle from the hole and gently proceed with the injection applying the least amount of pressure. Rely on the mark located on the hub indicating the orientation of the side port and the graduations every centimeter along the tube for determining the depth of the injection to favor a homogeneous treatment.
- 7/ To remove the SoftFil® micro-cannula from the syringe, hold the syringe and turn the mobile female luer lock counterclockwise firmly, then pull the micro-cannula backwards gently. Dispose of micro-cannulas and pre-hole needles in proper sharps container.

**SIDE EFFECTS**

Few complications as pain, redness, bruising, hematoma, swelling or oedema, may result from the injection technique itself. Furthermore, repeated injections can cause tissue scarring.

**VALIDITY / STORAGE CONDITIONS**

Do not use after the expiry date.

The product must be stored in its original packaging, in a dry place and in the shelter of light.

**SYMBOLS USED ON PACKAGING**

- Product reference   
 Use until   
 Batch number   
 Date of manufacture  
 Do not reuse   
 Do not resterilize   
 Do not use if package is opened or damaged  
 Sterile product. Method of sterilization – Ethylene oxide   
 Single sterile barrier system  
 See the Instructions for use leaflet   
 Unique Device Identifier   
 Dimension  
 Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician  
 Manufacturer   
 Keep away from moisture   
 Keep away from sunlight  
 Micro-cannula   
 Pre-hole needle   
 Latex free   
 Medical Device

**COLOUR CODES ON HUBS AND LABELS**

Ligh green = 14G    White = 16G    Pink = 18G    Yellow = 20G    Black = 22G    Dark blue = 23G  
 Medium purple = 24G    Orange = 25G    Brown = 26G    Middle grey = 27G    Yellow = 30G

All serious adverse events associated with use of this device should be reported to the manufacturer. All returned products must first be decontaminated, and accompanied by a decontamination certificate. If no certificate is provided, the claim will be rejected.

Patient Assistance Information: If you require immediate assistance, please contact your health care provider. For further assistance, please contact Prolenium Medical Technologies Inc. at 1-866-353-3017 (24 hours a day).

SoftFil® Precision, international patent-pending

Canada  
 Manufacturer / Fabricant:

United States  
 Manufactured For / Fabricado Para:

**Prolenium Medical Technologies Inc.**  
**138 Industrial Parkway North, Aurora**  
**Ontario L4G 4C3 CANADA**  
**T: 1.866 353 3015 | T: 1.905 508 1469**  
**www.prolenium.com**

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